

It's time to go universal.

Meet the universal series: Simplifying restorations from inventory to outcomes.

As a dentist, you are asked to do it all. You should have products that can do it all too – anterior and posterior, direct and indirect, self-cure and light cure¹, all classes, and all etching techniques.

3M is proud to offer a collection of universal products to simplify your day. These are all you need to perform virtually any restoration – so you can save time, reduce inventory, and streamline procedures with the confidence that comes from using world-class materials.

3M[™] RelyX[™] Universal Resin Cement works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M[™] Scotchbond[™] Universal Plus Adhesive.

3M[™] Scotchbond[™] Universal Plus Adhesive a true one-bottle system² with gold standard bond strength, suitable for all direct and indirect indications, all etching techniques, and all dental substrates. It also serves as a universal primer for all restoration materials.

3M[™] Filtek[™] Universal Restoratives simplify your single-shade restorations while providing natural-looking aesthetics with top-of-the-line polish retention and wear resistance.

Scotchbond Universal Plus Adhesive is dual cure compatible; Filtek composites are light cure only.
 No separate primer, silane, or dual-cure activator required.







Goodbye complexity. Introducing a truly universal resin cement.

3M[™] RelyX[™] Universal Resin Cement works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M[™] Scotchbond[™] Universal Plus Adhesive.

Together they eliminate the hassle and confusion of multiple resin cements, primers and adhesives, simplifying your direct and indirect restorative workflows. With fewer products to stock and replenish, you save space, time and money.

The Micro Mixing Tip is removed right after use enabling hygienic storage without a used mixing tip. The syringe is cleanly sealed by a unique valve mechanism.

Per application compared to currently available standard automix systems.
 Source: Field Evaluation EU/USA conducted by 3M.

The **Micro Mixing Tip** provides excellent mixing quality. Universal

Reit

Virtually all indirect indications with self-adhesive mode or as a twocomponent system with Scotchbond Universal Plus adhesive.

 Outstanding self-adhesion for maximum efficiency

Game-changing automix syringe with 80% less cement waste³

O Virtually no post-operative sensitivity⁴

Easy excess clean-up

Universal restoratives for every patient.

3M[™] Filtek[™] Universal Restoratives can offer versatility and simplicity for every class of restoration. They feature aesthetics and high polish retention for beautiful anterior restorations, while 3M's TRUE nanotechnology provides high strength and excellent wear resistance in posterior restorations.



3M[™] Filtek[™] Supreme XTE Universal Restorative

- Full palette of 36 shades and 4 opacities to match any patient's dentition
- Body shades enable a simple, single-shade technique
- Exceptional handling

Supreme XTE

- Wear resistance similar to enamel⁵
- Can be safely warmed⁶ and injection molded



3M[™] Filtek[™] Universal Restorative

- Simplified shade selection just 8 designer shades to cover most patients' aesthetic needs
- Universal opacity that blends naturally
- Extra White (XW) shade for patients who bleach
- Pink Opaquer to easily mask metal and stained dentition
- Can be safely warmed⁶ and injection molded

The original universal adhesive... Now better than ever!

First radiopaque universal adhesive

- Bonds and seals caries-affected dentine
- Gold standard adhesion to all dental substrates, including glass ceramics
- Universally suitable for all direct and indirect indications and all etching techniques
- Full dual- and self-cure compatibility without a separate dual cure activator

Virtually no postoperative sensitivity

3M Scotchb

Universal Plus Adhesive



Nearly a decade ago, 3M introduced the first truly universal dental adhesive. Now 3M[™] Scotchbond[™] Universal Plus Adhesive takes universal adhesives to the next level with new features like radiopacity for even greater control and predictability.

To simplify your workflow, Scotchbond Universal Plus Adhesive is suitable for all direct and indirect bonding indications and all etching techniques. It provides gold standard adhesion to all dental substrates, including glass ceramic. It is also a true one-bottle system – no separate primer or silane required – and it offers full dual and self-cure compatibility without a separate dual cure activator.

Also available in a unit dose for efficient hygiene management - reduces the risk of cross-contamination.

Ordering Information

3M[™] RelyX[™] Universal Resin Cement

Item No.	Product Information
Intro Kit	
56968	4 x 3M [™] RelyX [™] Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60 x 3M [™] RelyX [™] Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M [™] Scotchbond [™] Universal Plus Adhesive (5ml), 1 x 3M [™] Scotchbond [™] Universal Etchant (3ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides
Trial Kits	
56969	Shade TR
56970	Shade A1
	Each trial kit contains: 1 x 3M [™] RelyX [™] Universal Resin Cement syringe (3.4g), 15 x 3M [™] RelyX [™] Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M [™] Scotchbond [™] Universal Plus Adhesive (1.5ml), 1 x 3M [™] Scotchbond [™] Universal Etchant (3ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides
Refills	
56971	Shade TR
56972	Shade A1
56973	Shade A3O
56974	Shade WO
	Each refill contains: 1 x 3M [™] RelyX [™] Universal Resin Cement syringe (3.4g), 15 x 3M [™] RelyX [™] Universal Micro Mixing Tips, 5 Elongation Tips
Value Packs	
56977	Shade TR

56977 Shade TR 56978 Shade A1 Each value pack contains: 3 x 3M[™] RelyX[™] Universal Resin Cement syringes (3.4g each), 45 x 3M[™] RelyX[™] Universal Micro Mixing Tips, 15 Elongation Tips



3M.com.au/ScotchbondPlus



3M Australia Pty Limited Building A, 1 Rivett Rd North Ryde, NSW 2113 1300 363 454 www.3M.com.au

3M New Zealand Limited 94 Apollo Drive Rosedale, Auckland 0632 0800 808 182 www.3M.co.nz

3M[™] Scotchbond[™] Universal Plus Adhesive

Item No.	Product Information		
Vial			
41293	Intro Kit: 1 x 3M [™] Scotchbond [™] Universal Plus Adhesive vial (5ml), 50 disposable brushes, 1 x 3M [™] Scotchbond [™] Universal Etchant (3ml), 25 dispensing tips		
Vial Refill			
41294	1 x 5 ml vial		
41295	Refill Vial 3 x 5ml		
	3M™ Scotchbond™ Universal Plus Adhesive vial (5ml)		
Unit Dose			
41297	Intro Kit: 50 x L-Pop		
41298	100 x L-Pop		
41299	200 x L-Pop		
	3M [™] Scotchbond [™] Universal Plus Adhesive Unit Dose L-Pops, 1 x 3M [™] Scotchbond [™] Universal Etchant (3ml), 25 dispensing tips		

3M[™] Scotchbond[™] Universal Etchant Etching Gel

Item No.	Product Information	
41263	2 x 3M™ Scotchbond™ Universal Etchant (3ml each),	
	50 dispensing tips	
41274	3M [™] Scotchbond [™] Etchant Dispensing Tips	
	25 dispensing tips	

3M[™] Filtek[™] Universal Restorative

SYRINGE SHADE REFILLS - CAPSULES - include 20 - 0.2g capsules; SYRINGES - include 1 - 4g syringe

PINK OPAQUER (PO) REFILLS: CAPSULES – include 10 – 0.2g capsules; SYRINGES – include 1 – 4g syringe. All refills include a Technique Guide and Instructions for Use

Shade	Syringe	Capsule
Refills		
A1	6555A1	6550A1
A2	6555A2	6550A2
A3	6555A3	6550A3
A3.5	6555A3.5	6550A3.5
A4	6555A4	6550A4
B1	6555B1	6550B1
B2	6555B2	6550B2
D3	6555D3	6550D3
XW – Extra White	6555XW	6550XW
PO – Pink Opaquer	6555PO	

3M[™] Filtek[™] Supreme XTE Universal Restorative

SYRINGE SHADE REFILLS – includes 1 – 4g syringe; Instructions for Use CAPSULE SHADE REFILLS – Dentine, Enamel and Translucent Shade Refills include 10 – 0.2g capsules; Body Shade Refills include 20 – 0.2g capsules; Instructions for Use

Shade	Syringe	Capsule		
Body		·		
A1B	5914A1B	5915A1B		
A2B	5914A2B	5915A2B		
A3B	5914A3B	5915A3B		
A3.5B	5914A3.5B	5915A3.5B		
A4B	5914A4B	5915A4B		
A6B	5914A6B	5915A6B		
B1B	5914B1B	5915B1B		
B2B	5914B2B	5915B2B		
B3B	5914B3B	5915B3B		
B5B	5914B5B	5915B5B		
C1B	5914C1B	5915C1B		
C2B	5914C2B	5915C2B		
C3B	5914C3B	5915C3B		
D2B	5914D2B	5915D2B		
D3B	5914D3B	5915D3B		
WB	5914WB	5915WB		
XWB	5914XWB	5915XWB		
Dentine				
A1D	5914A1D	5915A1D		
A2D	5914A2D	5915A2D		
A3D	5914A3D	5915A3D		
A4D	5914A4D	5915A4D		
B3D	5914B3D	5915B3D		
C4D	5914C4D	5915C4D		
WD	5914WD	5915WD		
Enamel				
A1E	5914A1E	5915A1E		
A2E	5914A2E	5915A2E		
A3E	5914A3E	5915A3E		
B1E	5914B1E	5915B1E		
B2E	5914B2E	5915B2E		
D2E	5914D2E	5915D2E		
WE	5914WE	5915WE		
XWE	5914XWE	5915XWE		
Translucent				
Clear	5914CT	5915CT		
Blue	5914BT	5915BT		
Grey	5914GT	5915GT		
Amber	5914AT	5915AT		

3M, Filtek, RelyX and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. © 3M 2021. All rights reserved. Please recycle. AV011477567.